

ACCESSORIES PRICE LIST

ILLUSTRATIONS ARE NOT ACTUAL SIZE

Foam Telephone Pad...\$3.00

Bulb Air Blower\$6.00

Comply Ear Tipsea. \$4.00

Wax Picks\$2.50

Pouch\$4.00

Soft Wraps.....\$9.95/pkg. 10

BATTERY CHECKER/
HOLDER KEYCHAIN



\$3.95

TWIN-PACK BATTERY
CASE & CHAIN




\$1.50

FOR BEHIND-THE-EAR AIDS



\$5.00
Elbow

Please specify your Hearing Aid model



\$2.00
#13C Curved Tubing for Custom Earmolds



\$2.00/Ft.
#13T Straight Tubing



\$2.00
#TLCT Core & Tubing for Inserts

MUSHROOM
INSERTS



\$2.50
Each

TLM-SS ...X-Small
TLM-SSmall
TLM-M.....Medium
TLM-LLarge
TLM-LL....X-Large

DRI-AID KIT

\$8.00



Moisture is the enemy of hearing aids. Protect against corrosion and deterioration by removing the damaging moisture with Dri-Aid.

Inexpensive and very effective!
State choice of pouch or canister.

FOR BODY AIDS



Receivers:
Specify Model-\$30.00



\$2.00
#TLC Body Aid Connector



\$2.50
Mushroom Insert (See above)



CORDS

Straight Cord - \$20.00
Y-Cord - \$26.00

Please specify 2 or 3 prong and specify length and your aid model (18", 24", 30", or 36")

Is Your Hearing Aid
At Risk?



Ad-hearPlus
GERMINE GUARDS

The easiest way to keep your hearing aid clean

TROUBLE WITH WAX GETTING INTO YOUR AID? HERE'S THE ANSWER

Earwax getting inside the aid is by far the #1 cause for repairs of I-T-E hearing aids. Avoid costly repairs and downtime with your instruments. This proven remedy works similar to a "Band-aid" as it covers the receiver outlet hole with a special pass-through filter allowing the sound to pass clearly while keeping out damaging earwax or wetness. Each Ad•Hear lasts from 2-3 weeks. Available in different sizes... let us know for what model hearing aid you need them.

Ad•Hear

\$9.95/pkg. of 10

PLEASE SEND PAYMENT WITH ORDER Illinois residents - add 1% sales tax